

ROPOX

BETTER WAYS TO BETTER DAYS



**DOMINO
SLINGS**

For comfort and safety


DOMINO SLINGS

MADE OF GREY POLYESTER COMFORT MATERIAL

Domino Slings are produced using a specially manufactured polyester material, which unites comfort with an incredible strength. The material has an emery finish on the inside. This gives the Sling a cotton-like surface, which is both comfortable for the user and at the same time ensures that the Sling does not slide. The material is fire retardant.



DOMINO HYGIENE is incl. safety belt. Suitable for lifting persons from sitting to sitting position. Suitable for hygiene and bathroom purposes.



DOMINO ACTIVE. Supports the body up to shoulder height. Easy to use, as the leg and back sections are extra padded with special foam material. Safety belt as an accessory.



DOMINO COMBI is designed like the Active Sling, but has an extra built-in head support, which is hidden in a pocket in the back of the sling. Safety belt as an accessory.



DOMINO CLASSIC. Suitable for persons with little control over their head and body. Suitable for lifting from a lying position. Supports the whole body. Extra foam padding in leg, body and head section – for better comfort.



DOMINO AMPUTATION. The special design of the straps provides optimum support and balance. Safety belt as an accessory.



DOMINO THORAX. Only for stand-up hoist. Unique design and padding. Provides optimum support for stand-up lifting. Can also be supplied without padding under the arms.

DOLPHIN SLING – MADE OF WHITE POLYESTER NET FOR USE IN BATHROOMS

The water can easily pass through the material. Easy to clean.



DOLPHIN HYGIENE is incl. of safety belt. Suitable for hygiene and bathing purposes. Suitable for lifting persons from sitting to sitting position.



DOLPHIN CLASSIC. Suitable for persons with little control over their head and body. Suitable for lifting from a lying position. Supports the whole body.

DOMINO BLUE MADE OF SMOOTH BLUE POLYESTER MATERIAL (NYLON)


The low friction in the material makes it easier to place and remove the sling from the patient.



DOMINO BLUE HYGIENE is incl. of safety belt. Suitable for lifting persons from sitting to sitting position. Suitable for hygiene and bathroom purposes.



DOMINO BLUE ACTIVE. Supports the body up to the shoulder height. Safety belt as an accessory.



DOMINO BLUE CLASSIC. Suitable for persons with little control over their head and body. Suitable for lifting from a lying position. Supports the whole body.

ACCESSORIES



For the **ACTIVE, COMBI** and **AMPUTATION** Slings.
The use of Safety Belts gives the user a feeling of increased safety. They can easily be fixed to the slings.



Detachable 15cm and 30cm long strap loops are available. These can be fitted to the slings if you need increased possibility of individual adjustment and lifting position.

ORDERING INFORMATION

ART NO. DESCRIPTION

| DOMINO SLINGS – FOR STAND-UP HOIST | |
|------------------------------------|---|
| 25-43011 | Domino Thorax Sling with arm-padding. Small, max waist 70cm |
| 25-43012 | Domino Thorax Sling with arm-padding. Medium, max waist 100cm |
| 25-43013 | Domino Thorax Sling with arm-padding. Large, max waist 130cm |
| 25-43014 | Domino Thorax Sling with arm-padding. X Large, max waist 150cm |
| 25-43015 | Domino Thorax Sling with arm-padding. XX Large, max waist 170cm |
| 25-48011 | Domino Thorax Sling without padding. Small, max waist 70cm |
| 25-48012 | Domino Thorax Sling without padding. Medium, max waist 100cm |
| 25-48013 | Domino Thorax Sling without padding. Large, max waist 130cm |
| 25-48014 | Domino Thorax Sling without padding. X Large, max waist 150cm |

DOMINO SLINGS FOR BATH & TOILET

DOMINO HYGIENE

| | |
|----------|--|
| 25-48021 | Domino Hygiene Sling. Small, incl. safety belt, max waist 70cm |
| 25-48022 | Domino Hygiene Sling. Medium, incl. safety belt, max waist 110cm |
| 25-48023 | Domino Hygiene Sling. Large, incl. safety belt, max waist 125cm |
| 25-48024 | Domino Hygiene Sling. X Large, incl. safety belt, max waist 150cm |
| 25-48025 | Domino Hygiene Sling. XX Large, incl. safety belt, max waist 170cm |

DOMINO BLUE HYGIENE

| | |
|----------|--|
| 25-47011 | Domino Blue Hygiene Sling. Small, incl. safety belt, max waist 70cm |
| 25-47012 | Domino Blue Hygiene Sling. Medium, incl. safety belt, max waist 110cm |
| 25-47013 | Domino Blue Hygiene Sling. Large, incl. safety belt, max waist 125cm |
| 25-47014 | Domino Blue Hygiene Sling. X Large, incl. safety belt, max waist 150cm |

DOLPHIN HYGIENE

| | |
|----------|--|
| 25-42021 | Dolphin Hygiene Sling. Small, incl. safety belt, max waist 70cm |
| 25-42022 | Dolphin Hygiene Sling. Medium, incl. safety belt, max waist 110cm |
| 25-42023 | Dolphin Hygiene Sling. Large, incl. safety belt, max waist 125cm |
| 25-42024 | Dolphin Hygiene Sling. X Large, incl. safety belt, max waist 150cm |

DOMINO SLINGS WITH PADDED LEG- AND BACK SUPPORT

DOMINO ACTIVE

| | |
|----------|-----------------------------|
| 25-48041 | Domino Active Sling. Small |
| 25-48042 | Domino Active Sling. Medium |
| 25-48043 | Domino Active Sling. Large |

DOMINO BLUE ACTIVE

| | |
|----------|----------------------------------|
| 25-47041 | Domino Blue Active Sling. Small |
| 25-47042 | Domino Blue Active Sling. Medium |
| 25-47043 | Domino Blue Active Sling. Large |

ART NO. DESCRIPTION

| DOMINO SLINGS WITH FLEXIBLE HEAD SUPPORT AND PADDED LEG- AND BACK SUPPORT | |
|---|--|
| DOMINO COMBI | |
| 25-49031 | Domino Combi Sling. Small. Optional use of head support |
| 25-49032 | Domino Combi Sling. Medium. Optional use of head support |
| 25-49033 | Domino Combi Sling. Large. Optional use of head support |

DOMINO SLINGS WITH HEAD SUPPORT

DOMINO CLASSIC

| | |
|----------|---|
| 25-49011 | Domino Classic. Small. With padded leg supports |
| 25-49012 | Domino Classic. Medium. With padded leg supports |
| 25-49013 | Domino Classic. Large. With padded leg supports |
| 25-49014 | Domino Classic. X Large. With padded leg supports |

DOMINO BLUE CLASSIC

| | |
|----------|-----------------------------|
| 25-47841 | Domino Blue Classic. Small |
| 25-47842 | Domino Blue Classic. Medium |
| 25-47843 | Domino Blue Classic. Large |

DOLPHIN CLASSIC

| | |
|----------|--|
| 25-41011 | Dolphin Classic. Small. With padded leg supports |
| 25-41012 | Dolphin Classic. Medium. With padded leg supports |
| 25-41013 | Dolphin Classic. Large. With padded leg supports |
| 25-41014 | Dolphin Classic. X Large. With padded leg supports |

AMPUTATION SLINGS WITH PADDED LEG- AND BACK SUPPORT

DOMINO AMPUTATION

| | |
|----------|----------------------------------|
| 25-48031 | Domino amputations Sling. Small |
| 25-48032 | Domino amputations Sling. Medium |
| 25-48033 | Domino Amputations Sling. Large |

ACCESSORIES FOR DOMINO SLINGS

| | |
|----------|--|
| 25-49111 | Safety belt for Domino Active, Combi and Amputations Sling. Small |
| 25-49112 | Safety belt for Domino Active, Combi and Amputations Sling. Medium |
| 25-49113 | Safety belt for Domino Active, Combi and Amputations Sling. Large |
| 25-49121 | Safety belt for Domino Blue Active Sling. Small |
| 25-49122 | Safety belt Domino Blue Active Sling. Medium |
| 25-49123 | Safety belt for Domino Blue Active Sling. Large |
| 25-40115 | Strap loop 15cm. Can be used on all Domino Slings |
| 25-40130 | Strap loop 30cm. Can be used on all Domino Slings |

ROPOX

BETTER WAYS TO BETTER DAYS

Ropox A/S

Ringstedgade 221

DK-4700 Naestved

Tel. +45 5575 0500

Fax +45 5575 0550

info@ropox.dk

www.ropox.com